

## i-DIGITS™ QUANTUM

Design, dexterity and  
intelligent motion

The first individually designed, multi-articulating partial hand prosthesis that can be controlled with simple gestures.

- Suitable for individuals with partial hand absence of 1-5 digits at or proximal to the metacarpophalangeal joint.
- Digits move independently and work in conjunction with any remaining fingers.
- Gesture Control powered by i-Mo™ technology uses simple gestures to change grips
- Bluetooth enabled Grip Chips™ activate pre-set grips in proximity to i-Digits
- Speed boost:
  - Adjustable open/close speed for digits



i-Mo™ Intelligent Motion  
GESTURE CONTROL



Triggers  
MUSCLE CONTROL




Quick Grips™  
APP CONTROL





Grip Chips™  
PROXIMITY CONTROL



i-Digits Quantum is appropriate for partial hand absence where the level of loss or deficiency is distal to the wrist and proximal of the metacarpophalangeal joint. Anywhere from one to five digits can be replaced.





### Digit mounting options



XS mounting plate: straight digit position

M mounting plate: splayed digit position

### Digit size options:



### i-DIGITS QUANTUM

Part #	Description
TBX50008	PK1 i-Digit Quantum 1 Digit
TBX50016	PK2 i-Digit Quantum 2 Digit
TBX50024	PK3 i-Digit Quantum 3 Digit
TBX50032	PK4 i-Digit Quantum 4 Digit
TBX50040	PK5 i-Digit Quantum 5 Digit

### PACKAGE CONTENTS\*

Electrode kit	yes
Grip Chips™	yes
Battery and charger	yes
Car charger	yes
Socket components/mounting plates	yes
iPod Touch	yes
i-Limb skin coverings	3 coverings per finger

\*2 years of touch care warranty

FOLLOW ÖSSUR ON



USA (800) 233-6263  
WWW.OSSUR.COM

CANADA (800) 663-5982  
WWW.OSSUR.CA

