

REBOUND® DUAL ST

Evolving from the functional support and versatility of the Rebound® DUAL, Össur offers a knee ligament brace for everyone. The Rebound DUAL ST, a sleek step-through design built with the strength of aircraft aluminum. Leveraging a light-weight profile and patented hi-tech ActiveGrip™ liners, the Rebound DUAL ST is designed for exceptional comfort and fitting with ease.



KNEE


INDICATIONS FOR USE

For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, and combined instabilities

FEATURES AND BENEFITS

- Step-Through Design: advantage in comfort & fitting with ease
- Sleek, Strong Frame: offering the strength of aircraft aluminum with a light-weight design, the 14" (35 cm) frame length keeps a low-profile
- Anti-Gravity: Activegrip™ liners with the added advantage of Össur's anti-migration strap (AMS) system
- Easy for Everyone: simple strapping system and field-serviceable technology for adjustments on the go
- Versatile, comfortable, and reliable support for improved ligament stability
- Suitable for low-to-high impact activities as well as sliding sports
- Optional PCL kit accessory for support of the PCL
- Available in pre-fabricated sizes XS-3XL and also custom-made

REQUIRED MEASUREMENTS



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Knee caliper		Thigh circumference 15 cm (6") above mid patella	
	cm	inch	cm	inch
XS	7.6 - 8.9	3 - 3½	38.1 - 41.9	15 - 16½
S	8.9 - 10.2	3½ - 4	41.9 - 47	16½ - 18½
M	10.2 - 11.4	4 - 4½	47 - 52.1	18½ - 20½
L	11.4 - 12.7	4½ - 5	52.1 - 57.2	20½ - 22½
XL	12.7 - 14	5 - 5½	57.2 - 62.2	22½ - 24½
XXL	14 - 15.2	5½ - 6	62.2 - 69.9	24½ - 27½
3XL	15.2 - 16.5	5½ - 6½	69.9 - 76	27½ - 30

REBOUND® DUAL ST

ORDERING INFORMATION

PRE-FABRICATED

Size	Left	Right
XS	B-242620001	B-242520001
S	B-242620002	B-242520002
M	B-242620003	B-242520003
L	B-242620004	B-242520004
XL	B-242620005	B-242520005
XXL	B-242620006	B-242520006
3XL	B-242620007	B-242520007

ACCESSORIES

PCL STRAP KIT

Size	Part#
XS-L	B-701090065
XL-3XL	B-701090066

CALF LINER

Size	Part#
XS	B-701090074
S	B-701090075
M	B-701090076
L	B-701090077
XL	B-701090078
XXL	B-701090079
3XL	B-701090080

THIGH LINER

Size	Part#
XS	B-701090067
S	B-701090068
M	B-701090069
L	B-701090070
XL	B-701090071
XXL	B-701090072
3XL	B-701090073

FLEXION STOP KIT

Description	Part#
Flexion stop kit	B-701090004

STRAP KIT

Size	Part#
XS-L	B-701090005
XL-3XL	B-701090006

OTHER ACCESSORIES

Description	Part#
Extension stop kit	B-701090008
Condyle pad kit	B-701090009
Field servicable D-ring kit	B-701090059

REBOUND® DUAL ST CUSTOM

Evolving from the functional support and versatility of the Rebound® DUAL, Össur offers a knee ligament brace for everyone. The Rebound DUAL ST, a sleek step-through design built with the strength of aircraft aluminum. Leveraging a light-weight profile and patented hi-tech ActiveGrip™ liners, the Rebound DUAL ST is designed for exceptional comfort and fitting with ease.



25 Color Options*

Customize to match your team colors



*Custom models only. Colors are approximate

KNEE


INDICATIONS FOR USE

For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities


FEATURES AND BENEFITS

- Step-Through Design: advantage in comfort & fitting with ease
- Sleek, Strong Frame: offering the strength of aircraft aluminum with a light-weight design, the 14" (35 cm) frame length keeps a low-profile
- Anti-Gravity: Activegrip™ liners with the added advantage of Össur's anti-migration strap (AMS) system
- Easy for Everyone: simple strapping system and field-serviceable technology for adjustments on the go
- Versatile, comfortable, and reliable support for improved ligament stability
- Suitable for low-to-high impact activities as well as sliding sports
- Optional PCL kit accessory for support of the PCL
- Custom-made (also available in pre-fabricated sizes XS-3XL)


REQUIRED MEASUREMENTS




Measure calf circumference 15 cm (6") below mid-patella. Leg in extension.



Measure calf circumference 7.5 cm (3") below mid-patella. Leg in extension.



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.



Measure thigh circumference 7.5 cm (3") above mid-patella. Leg in extension.

REBOUND® DUAL ST CUSTOM

ORDERING INFORMATION

CUSTOM

Left	Right
B-155620002	B-155520001

ACCESSORIES

FLEXION STOP KIT

Description	Part#
Flexion stop kit	B-701090004

PCL STRAP KIT

Size	Part#
XS-L	B-701090065
XL-3XL	B-701090066

SPORT UNDERSLEEVE

Size	Part#
S	20682
M	20683
L	20684
XL	20685

NEOPRENE UNDERSLEEVE

Size	Part#
S	20251
M	20252
L	20253
XL	20254
XXL	20255
3XL	20256
4XL	20257

NEOPRENE OVERSLEEVE

Size	Part#
S	20051
M	20052
L	20053
XL	20054
XXL	20055
3XL	20056
4XL	20057

PADDED SPORTS OVERSLEEVE

Size	Part#
S	20448
L	20449

CALF LINER

Size	Part#
XS	B-701090074
S	B-701090075
M	B-701090076
L	B-701090077
XL	B-701090078
XXL	B-701090079
3XL	B-701090080

THIGH LINER

Size	Part#
XS	B-701090067
S	B-701090068
M	B-701090069
L	B-701090070
XL	B-701090071
XXL	B-701090072
3XL	B-701090073

STRAP KIT

Size	Part#
XS-L	B-701090005
XL-3XL	B-701090006

OTHER ACCESSORIES

Description	Part#
Extension stop kit	B-701090008
Condyle pad kit	B-701090009
Field servicable D-ring kit	B-701090059