



LIVINGSKIN[™]

Passive functional prostheses made with high definition silicone.

- Custom sculpting & painting ensure an individualized fit
- Mimics a realistic appearance by matching skin tone and adding freckles, hair, scars and tattoos.
- The functional capabilities of the Livingskin prosthesis can assist in Activities of Daily Living such as pushing, pulling, stabilizing, supporting and typing.
- Skin color and other features are added below the outer clear silicone layer to prevent fading.
- Adjustable armatures may be fitted to allow repositioning of each finger for specific tasks.
- Silicone provides a versatile, durable and virtually stain-proof prosthetic option.

NEW & IMPROVED

- Ask us about our brand new pediatric offering!


VARIOUS OPTIONS AVAILABLE:



ACCEPTANCE GUARANTEE

Össur is committed to deliver the highest quality Livingskin products. Our Acceptance Guarantee ensures that the product is delivered to meet your patient's needs.

Our revolutionary iPad-based color capture system enables you to accurately match skin tones and select features such as hair, nails, veins, or freckles. With an intuitive design, the easy on screen step-by-step instruction guides you through the entire process resulting in a high quality prosthesis.



Livingskin Fingers

| Part # | Part # (Pediatric) | Description |
|--------|--------------------|-------------|
| 066026 | 066704P | 1 Finger |
| 066029 | 066705P | 2 Finger |
| 066039 | 066706P | 3 Finger |
| 066040 | 066707P | 4 Finger |
| 066134 | 066708P | 5 Finger |
| 066135 | 066714P | Finger Band |

Livingskin Hands

| Part # | Part # (Pediatric) | Description |
|--------|--------------------|--------------|
| 066014 | 066709P | Partial Hand |
| 066015 | 066710P | Full Hand |

Livingskin Arms

| Part # | Part # (Pediatric) | Description |
|--------|--------------------|----------------------------------|
| 066703 | 066711P | Wrist Disarticulation Passive |
| 066016 | 066712P | Transradial Passive to Elbow |
| 066017 | 066713P | Transhumeral Passive to Axilla |
| 066454 | N/A | Shoulder Disarticulation Passive |

i-Limb Skin Match Gloves*

| Part # | Description |
|--------|---|
| 066006 | i-Limb Skin Match Hand Cover |
| 066448 | I-Limb Skin Match Forearm Sleeve |
| 066449 | I-Limb Skin Match Hand Covering including Custom Forearm Sleeve |

* PASSIVE UPPER LIMB PROSTHETICS LIMITED WARRANTY - Livingskin Passive Device - 12 Months

Livingskin Toes

| Part # | Part # (Pediatric) | Description |
|----------|--------------------|-------------|
| 066130A | 066715P | 1 Toe |
| 066131A | 066716P | 2 Toe |
| 066132A | 066717P | 3 Toe |
| 066133A | 066718P | 4 Toe |
| TBX50404 | 066719P | 5 Toe |
| TBX50410 | 066714P | Toe Band |

Livingskin Feet

| Part # | Part # (Pediatric) | Description |
|----------|--------------------|--------------|
| TBX50374 | 066720P | Partial Foot |
| 066136A | 066721P | Full Foot |

Livingskin Accessories

| Part # | Description |
|----------|--|
| TBX50160 | SeasonGuard |
| 066443 | Dermahair - Finger |
| 066065 | Dermahair - Hand |
| 066442 | Dermahair - Arm |
| 092039 | Extra Silicone Material - per 1 lb Tube |
| TBX50181 | Orthotic Device (including Thumb Stabilizer) |
| TBX50135 | Pronation/Supination Wrist Component |
| TBX50070 | Flexion/Extension Wrist Component |
| TBX50166 | Slipcoat |

FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM

CANADA (800) 663-5982
WWW.OSSUR.CA

